

CyberCryo Cold Trap

- Cold Trap capacity 2 or 4 or 6 liters
- Condensing capacity from 750 to 450 ml
- Vacuum gauge included
- Max. temperature -60°C
- Fuzzy free
- CFC free
- Noiseless



Accessories PVD Closure plate with vacuum in and outlet GE Glass insert trap for condensing isotop material

MC 2 L —60°C

Danfoss 1/6 PS Compressor

CFC free system

Freon R507-x 70gr

Power 220 / 240V

Amp / Watt 1.8 / 400

Condenser surface 1460cm²

Condenser made in stainless steel

Outside dim. 290 x 535 x 265

Weight 20 kg

MC 4 L —60°C

Danfoss 1/5 PS

CFC free system

Freon R507-x 90 gr

Power 220 / 240V or 110V

Amp / Watt 2 / 450

Condenser surface 1460cm²

Condenser made in stainless steel

outside dim. 400 x 685 x 255

Weight 32 kg

MC 6 L —60°C

Danfoss 1/5 PS

CFC free system

Freon R507-x 100gr

Power 220 / 240V or 110V

Amp / Watt 2 / 450

Condenser surface 1460cm²

Condenser made in stainless steel

Outside dim. 400 x 700 x 335

Weight 35 kg

Optional Accessory

CyberCryo Freeze Dryer

Table Top Equipments

- Condensing Vessel Capacity 2 or 4 or 6 liters
- from 750 up to 4500 ml condensing capacity
- Basic Freeze Dryer consisting of defrost heater and vacuum gauge
- Max. temperature -60°C
- FCKW free
- Fuzzy logic



MC 2L With Disc Manifold



MC 4L With Drum Dryer

- Transparent Freeze Drying drum
- Different size for 2; 4 and 6 liters Freeze drying condenser
- Variable drying platforms with heated or non heated surface
- Power supply for heating platform 24 Volt
- Special design on request
- Accessories Optional



MC 2L -60°C



MC 4L -60°C



MC 6L -60°C

Cyberlab Corporation
22, Salo Terrace,
Millbury, MA01527
USA
Tel.: 001 508 2273494
Fax: 001 806 2091960

Cyberlab (Asia-Pacific) Pte. Ltd.
9, Amber Road,
#16-02, Parkway Mansion,
Singapore 439954
Tel.: 0065 63397696
Fax: 0065 63370048

Cyberlab Gulf Ltd.
P. O. Box # 261597,
Jebel Ali FZ, Dubai,
U.A.E.
Tel.: 00971 6 5738582
Fax: 00971 6 5738584

Email : info@cyberlabusa.com
URL : www.cyberlabusa.com