


# Rotary Evaporators

CYBERLAB™

*CYBERLAB ALWAYS THINKS FOR CUSTOMER SATISFACTION*



## Rotary Evaporator **CR-2000**



- 6 types of flasks ranging from 50ml to 3,000ml
- Cyberlab is equipped with speed control and R.P.M. value is displayed in digital.
- Digital Water bath covers ambient temperature up to 180°C
- If water level at heating bath goes down, power will be cut off automatically by sensor.
- Speed Control: 20~280 R.P.M.
- Slant condenser
- Standard distillation
- Continuous supply device for large capacity
- Central connection
- Dimensions (WxDxH): 320x730x520 mm
- Power Source: AC 220V 50/60 Hz 2.5A
- Heating Bath Power Wattage: 1,000W 4.5A
- Weight: 16kg

## Rotary Evaporator **CR-2001**



- 6 types of flasks ranging from 50ml to 3,000ml
- Cyberlab is equipped with speed control and R.P.M. value is displayed in digital.
- Digital Water bath covers ambient temperature up to 180°C
- If water level at heating bath goes down, power will be cut off automatically by sensor.
- Speed Control: 20~280 R.P.M.
- Vertical condenser
- Standard distillation
- Continuous supply device for large capacity
- Central connection
- Dimensions (WxDxH): 320x580x720 mm
- Power Source: AC 220V 50/60 Hz 2.5A
- Heating Bath Power Wattage: 1,000W 4.5A
- Weight: 17kg



### **Cyberlab Corporation**

22, Salo Terrace,  
Millbury, MA01527  
USA  
Tel.: 001 508 2273494  
Fax: 001 806 2091960

### **Cyberlab (Asia-Pacific) Pte. Ltd.**

9, Amber Road,  
#16-02, Parkway Mansion,  
Singapore 439954  
Tel.: 0065 63397696  
Fax: 0065 63370048

### **Cyberlab Gulf Ltd.**

P. O. Box # 261597,  
Jebel Ali FZ, Dubai,  
U.A.E.  
Tel.: 00971 6 5738582  
Fax: 00971 6 5738584

Email : [info@cyberlabusa.com](mailto:info@cyberlabusa.com)  
URL : [www.cyberlabusa.com](http://www.cyberlabusa.com)