

Ultrasonic Cleaner

Bench-top Ultrasonic Cleaner

Ultrasonic cleaning is faster, more thorough and safer than any other method-scrubbing, soaking, steam and so on. Ultrasonic sound waves moving through a cleaning solution create an effect called cavitation, the rapid formation and implosion of microscopic bubbles. Heat enhances the process. The deep cleaning action of ultrasonics removes even the most stubborn dirt and contaminants.



Features

- Powerful and reliable cleaning
- Built - in microprocess controller
- Furnishing all information - ultrasonic power, set time, elapsed time, set temperature, current temperature, heater operation status and so on for optimum conditions / operation status.
- Selectable use with ultrasonic strength controller(3steps-L.M.H)
- Built-in a digital thermometer and timer to set a optimum condition for cleaning
- Easy to read operation status on digital display(400 series)
- Built - in drain valve unit - Three model case on engineering plastic molds
- Using a silicon rubber heater for durability and stability

Model	405	410	420
Bath Capacity (ℓ)	5.7	10.0	20.7

CyberBath-CB2150 Ultrasonic Cleaner



Internal tank size	295mm x 150mm x 100mm
External dimension	320mm x 175mm x 240mm
Capacity	4.5 liters
Material	Internal Tank: S. S.#304 External Body: Aluminium 5052 Cover: S.S.#304 Basket: S.S. #304
Ultrasonic output	200W / 46KHZ
Heater output	150W
Electronic timer	1-30 minutes (adjustable)
Elec. temp. con.	Ambient temp.-60°C (Adjustable)
Power switch	Soft touch button ON/OFF switch
Power source	110V, 220V, 240V as order request
Standard accessories	S.S. cover x 1Pc. & S.S. basket x 1Pc

CyberBath-CB2080 Ultrasonic Cleaner



Internal tank size	150mm x 135mm x 100mm
External dimension	175mm x 160mm x 210mm
Capacity	2 liters
Material	Internal Tank: S. S.#304 External Body: Aluminium 5052 Cover: S.S.#304
Ultrasonic output	80W / 46KHZ
Electronic timer	1-30 minutes (adjustable)
Power switch	Soft touch button ON/OFF switch
Power source	110V, 220V, 240V as order request
Standard accessories	S.S. cover x 1Pc.

CyberBath-CB2050 Ultrasonic Cleaner



Internal tank size	168mm x 89mm x 55mm
External dimension	193mm x 116mm x 118mm
Capacity	0.6 liters
Material	Internal Tank: S. S.#304 External Body: ABS Plastic (U/L fire resistant) Cover: PS plastic Basket: ABS plastic (U/L fire resistant)
Ultrasonic output	50W / 46KHZ
Electronic timer	Un-adjustable
Power switch	Soft touch button ON/OFF switch
Power source	110V, 220V, 240V as order request
Standard accessories	Plastic cover x 1Pc. & plastic basket x 1Pc.
Safety approval	CE marking