

Portable Electronic Balance

US-Series

Precision without compromise in the field



Features:

- Designed for accuracy
- High Precision
- Fast Response
- Large Pan Size
- LCD Display with Backlight
- AC & Battery Operations
- In-use Protective Cover
- Easy to operate
- Wide range of capacity, 300g - 10kg
- High Accuracy, 10mg/100mg
- Screw mounted protection for pan support shaft
- Economical & affordable Price
- 1 Year Warrantee

Applications:

- Gems & Jewellery
- Laboratory
- Banks
- Education
- Food Processing
- Retail

Specifications

Model	High resolution type			Standard resolution type		
	US-300g	US-600g	US-1000g	US-3000g	US-6000g	US-10000g
Capacity	300 g	600 g	1000 g	3000 kg	6000 g	10000 g
Minimum display	0.01 g	0.01 g	0.01 g	0.1 g	0.1 g	0.1 g
Taring range	0 g to capacity					
Ambient temperature	5 - 38 degree centigrade					
Unit conversion	g, ct, oz, dwt					
Data output	DATA I/O IFB interface connector					
Pan size	135mm dia		160 x 180mm			
Dimensions (approx.)	240W x 180D x 55H mm					
Power source	6 AA batteries, or AC adapter 6V					
Standard accessories	Main Balance, Weighing Pan, AC adapter, In-use protective cover					

• Secondary / Remote LED Display Optional

• Acrylic wind Shield Optional

Your Local Representative

Cyberlab Corporation

22, Salo Terrace,
Millbury, MA01527
USA
Tel.: 001 508 2273494
Fax: 001 806 2091960

Cyberlab (Asia-Pacific) Pte. Ltd.

9, Amber Road,
#16-02, Parkway Mansion,
Singapore 439954
Tel.: 0065 63397696
Fax: 0065 63370048

Cyberlab Gulf Ltd.

P. O. Box # 261597,
Jebel Ali FZ, Dubai,
U.A.E.
Tel.: 00971 6 5738582
Fax: 00971 6 5738584

Email : info@cyberlabusa.com
URL : www.cyberlabusa.com