

Chromatography Columns & Supplies

CyberClouser - Autosampler Vials, Caps, Septas & Inserts

Screw Top Vials, Caps & Septa

Combine the excellent autosampler handling of a crimp cap profile with the ease of a screw cap. The screw caps and precision fit septa give a secure seal with the microvolume inserts.

Universal Autosampler Vials (1000/pk)

- ST-09-2067 Screw Top Vial, 1.5mL, 12x32mm, 9mm, Clear
- ST-09-2068 Screw Top Vial, 1.5mL, 12x32mm, 9mm, Clear with Label
- ST-09-2069 Screw Top Vial, 1.5mL, 12x32mm, 9mm, Amber with Label
- ST-PP-1205 Screw Top Vial, 1.5mL, 32x11.6mm, 9mm, PP transparent with filling lines

- ST-08-4428 Screw Top Vial, 1.5mL, 12x32mm, 8mm, Clear
- ST-08-4429 Screw Top Vial, 1.5mL, 12x32mm, 8mm, Amber
(Micro- Inserts also available for 8mm Vials)

- CT-11-4491 Crimp Top Vial, 1.5mL, 12x32mm, 11mm, Clear
- CT-11-4492 Crimp Top Vial, 1.5mL, 12x32mm, 11mm, Clear with Label
- CT-11-4493 Crimp Top Vial, 1.5mL, 12x32mm, 11mm, Amber with Label

- ST-13-0222 Screw Top Vial, 4ml , 45x14.75mm,13mm,Clear Glass
- ST-13-1335 Screw Top Vial, 4ml , 45x14.75mm,13mm,Clear Glass with Label
- ST-13-0280 Screw Top Vial, 4ml , 45x14.75mm,13mm,Amber Glass
- ST-13-1336 Screw Top Vial, 4ml , 45x14.75mm,13mm,Amber Glass with Label
(0.3ml Micro-Insert with Spring also available for these Vials)

Micro Inserts

- MI-09-0357 0.1ml Micro-Insert,31*6mm,clear glass,15mm top(for 9mm / 11mm vial)
- MI-09-0865 0.1ml Micro-Insert,28*6mm,clear glass,with assembled plastic Spring(for 9mm/11mm vial)
- MI-09-0866 0.2ml Micro-Insert,31*6mm,clear glass,flat bottom (for 9mm/11mm vial)

Shell Vials with PE Plug

- SV-01-0513 1ml Clear Glass Shell Vial with PE Plug Cap,40x8mm
- SV-01-0595 1ml Amber Glass Shell Vial with PE Plug Cap,40x8mm
- SV-02-0544 2ml Clear Glass Shell Vial with PE Plug Cap,31.5x11.6mm
- SV-02-0545 2ml Amber Glass Shell Vial with PE Plug Cap,31.5x11.6mm
- SV-04-0548 4ml Clear Glass Shell Vial with PE Plug Cap,44.6x14.65mm
- SV-04-0562 4ml Amber Glass Shell Vial with PE Plug Cap,44.6x14.65mm
(Micro-Inserts for the above also available)

Universal Caps for 1.5ml Vials (1000/pk)

- SC-09-0726 PP Red Colour Open Top Screw Cap, 6mm hole w/o septa
- SC-09-0727 PP Green Colour Open Top Screw Cap, 6mm hole w/o septa
- SC-09-0728 PP Blue Colour Open Top Screw Cap, 6mm hole w/o septa

- SC-08-0728 PP Blue Colour Open Top Screw Cap, 5.5mm hole w/o septa
- SC-08-4438 PP Black Colour Open Top Screw Cap,5.5mm hole w/o septa

- CC-11-4498 11mm Silver Aluminum Crimp Caps with PTFE/Natural Rubber Septa
- CC-11-4500 11mm Silver Aluminum Crimp Caps with PTFE/Silicone Septa

- SC-13-0166 PP Screw Cap black,8.5mm Centre hole.

Universal Septa for Screw Caps (1000/pk)

- SP-09-0731 Ivory PTFE/Red Rubber Septa 9mm
- SP-09-0730 Red PTFE/White Silicone Septa 9mm
- SP-09-2074 Blue PTFE/White Silicone Septa with Slit for Waters

- SP-08-4437 Red PTFE/White Silicone Septa 8mm
- SP-13-0730 Silicone White /PTFE red Septa, 13mm

