

CYBERLAB™

Cyber-pH

pH meter

Model:
Cyber pH-14L



| Catalogue number | Model | Description |
|------------------|-------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------|
| CYB-pH14L | Cyber pH14L | Meter + tri-combination pH/ATC electrode + standard buffer solution (pH 4.00, 7.00, 10.01 x 50ml for each bottle) + standard electrode holder + AC/DC adaptor |

Features:

- Smart single chip design
- Apheliotropic LCD, display both pH and temperature value at the same time.
- Automatic calibration: 5 point (pH 4.00, 6.86, 7.00, 9.18 and 10.01) can be updated.
- Display slope percent automatically after calibrating
- °C and °F can be switched automatically.
- 25 datum can be stored or picked and can record max. and min. value.



Specification:

- pH range: 0~14.00 pH
- Resolution: 0.1/0.01 pH
- Accuracy: ±0.01 pH ±1 display digit
- Temperature compensation: 0~100 Auto/Manual
- Data connector: RS232 output
- Operation function: Automatic calibration, data storage, record max. and min. value
- Dimension and weight: 160 x 190 x 70mm/ 750g
- pH electrode: 201B-F plastic pH tri-combine electrode (include ATC)
- Temp. sensor (optional)
- 601 stirring frame (optional)
- pH buffer solution: A bottle of pH 4.00, 7.00, 10.01 respectively
- pH buffer reagent: A bag of pH 4.00, 7.00, 10.01 respectively



Cyberlab Corporation

22, Salo Terrace,
Millbury, MA01527
USA
Tel.: 001 508 2273494
Fax: 001 806 2091960

Cyberlab (Asia-Pacific) Pte. Ltd.

9, Amber Road,
#16-02, Parkway Mansion,
Singapore 439954
Tel.: 0065 63397696
Fax: 0065 63370048

Cyberlab Gulf Ltd.

P. O. Box # 261597,
Jebel Ali FZ, Dubai,
U.A.E.
Tel.: 00971 6 5738582
Fax: 00971 6 5738584

Email : info@cyberlabusa.com
URL : www.cyberlabusa.com